

Ai-surgery

The Latest Generation Piezo Bone Surgery

STRONG POWER

Strong power, no handpiece overheating, thanks to:

- New enhanced handpiece
- Advanced sensing technology
- Smart control system

Smart control system

- Automatic identification of different bone quality and strength to perform power adjustment, achieving the optimal output
- Easy and effort-saving for doctors
- Less noise at work



INTELLIGENT AND POWERFUL

Product with the latest technology from WOODPECKER, providing dentists with a superbly reliable and safe ultrasonic surgical system

- Sinus lift & augmentation surgery
- Bone surgery
- Bone ridge reduction surgery
- Periodontal surgery
- Implant surgery
- Endo surgery
- Maxillomandibular advancement surgery
- Cosmetic dental surgery

EXTREMELY RELIABLE

Convenient operation

Fully upgraded Woodpecker fourth generation of Piezo Bone Surgery, easier to operate



Mode selection

Power varies as per different bone quality



One-touch enhancement

Exclusively featured one-touch power enhancement



Water flow adjustment

Different water flow adjustment levels to choose from



Easy rinse

After each use, rinse the pipeline via control of foot pedal



Multi-function foot pedal

- Control of multiple functions
- One-button rinse mode



Sterilizable water pipe

- External handpiece for easy disassembly and maintenance
- Autoclavable

Ultra-quiet peristaltic pump

- Continuous and precise flow control; direct water flow to the tip end; efficient cooling
- More silent and stable






SURGERY TIPS



SC Series






SE Series




SL Series



SP Series



SS Series



SURGERY TIPS

SS1TR




SS11L



SS18



SS21




SS30



SS31



SS32



SS33



SS34




Basic Kits


SI1




SC1



SS1L



SS1R




SS1




SS2



SS3




SS4



SS5



SL1



SL2



SL3



SL4



SL5

