



SKR054-ET

EMERGENCY TROLLEY

Register code: SKM-A



SKR054-ET

FEATURES AND HIGHLIGHTS


Technical parameters

- External size (LxWxH) 750x520x960 mm
- Castors 100 mm



One-piece countertop

The one-piece countertop provides more comprehensive protection; the fully guarded handle makes it more comfortable to hold me; the raised grip shape makes the operation of the side-drawing board more comfortable.




Decorative gripper

Outer beveled decorative grip, better operating comfort and more three-dimensional shape.



Heavy duty construction

Wheels assembled on the same vertical line as four columns. More stable, even the drawer open inside full of medicine/tools.



Drawer

Drawer in size big & deeper!



Power supply

The power socket is movable and detachable, easy to use, located above the CPR board to convenient the use of oxygen cylinders




Inserted oxygen bottle holder

Inserted retractable oxygen bottle holder, belts fixed on trolley main body, space saving



Infusion stand

IV pole inserted into main body column.



Drawer holder

Card insert from downside, easy keep clean.



Defibrillator

The plate surface is made of virgin ABS material. The position of the I,V pole and the defibrillator can be exchanged for easy use.



Sorting trash can

Green+yellow sorted trash can.





CONSTANT
INNOVATION



Design



Research



Manufacturing



Maintenance



Service

Distribuidor

Beracah
MEDICA

✉ ventas@beracahmedica.mx

☎ 662 15 10 979